

Self-contained mobile spectrometer
TRISTAN 5

efficient photonic solutions

深圳市新世联科技有限公司


TRISTAN 5



Self-contained mobile spectrometer

The TRISTAN 5 is a compact and totally self-contained spectrometer with an integrated display. It is controlled via a color touch-screen display with a user interface similar to Windows or a PDA. Spectroscopy has never been so mobile and easy-to-handle before.

TRISTAN 5 comes as a complete package solution, containing everything you need for mobile spectroscopy: From the fully-equipped optical cell with entry slit to the calibration at our facility and the software package for professional spectroscopy.



- 14x2048 Pixel back-thinned CCD sensor
- Rechargeable battery
- Touch-screen color display

Latest sensor technology

The TRISTAN 5 contains a back-thinned CCD sensor with 14 x 2048 pixels. Compared to regular CCD arrays, back-thinned sensors offer much higher sensitivity in UV and NIR region without need for additional filters and coatings. The four standard configurations of TRISTAN 5 set new standards for off-the-shelf spectrometers with regard to sensitivity, dynamics and noise and make it a professional measuring device in the field as well as for the lab.


Onboard data processing

TRISTAN 5 is ideally suited for spectral analysis in the field. TRISTAN 5 can operate without external power supply for about four hours with the integrated battery pack. Measurement results are displayed right away. The onboard software controls the spectrometer with all its functions including the ready-to-use data preparation. Processing data and storage is performed onboard within the spectrometer. The internal memory of 4 MB can store up to 100 spectra, and an additional memory card can increase storage to 128 MB.

PC operation by TRIWin 3 Pro

TRISTAN 5 is controlled via a color touch-screen display with a user interface similar to Windows or a PDA. It comes with an USB interface and a RS232 serial port and can be upgraded with an Ethernet port. The included TRIWin 3 Pro software provides full control of your spectrometer from any Windows-based PC. TRISTAN 5 is the perfect spectrometer in the field as well as for the lab.

When connected to a PC and operated by TRIWin the TRISTAN 5 becomes a professional desktop spectrometer without any limitations. It is operated through an intuitive Graphical User Interface which offers an unrivaled abundance of tools for working with collected data, from archiving and documenting to performing mathematical functions for data analysis and curve smoothing. Output can be exported directly to Excel or saved in graphical format, which makes documentation very easy. Separate laboratory software is not necessary.



TRIWin Software

深圳市新世联科技有限公司

Options



Back-thinned CCD detector

Customized optical engine

The performance of any spectrometer is determined by the combination of wavelength range, resolution and sensitivity. Off-the-shelf spectrometers are targeted to the most frequently requested configurations. To meet your specific requirements the optical engine of all our TRISTAN spectrometers with back-thinned CCD arrays can be configured within a wide variety of specifications.

Please ask your local sales representative for further details.

Ethernet Port

One main feature of TRISTAN spectrometers is their stand-alone capability. With their onboard processing unit they can work on their own. This feature makes them the perfect device for remote measurements via LAN. With the Ethernet option TRISTAN can communicate with every PC in your LAN directly. The Ethernet port replaces the USB port in this design.

2nd Order Filter

Every dispersive element produces a signal in the desired direction, and it also produces a small reflection of this signal in the direction of higher wavelengths. This effect is called second order signal. For spectrometers with high bandwidth it is possible that such a second order signal might appear within the measured wavelength range. m-u-t offers second order filters to reduce this effect mounted directly on the surface of the detector.

Accessories

m-u-t offers a wide range of accessories for your spectroscopic measurement setup. These accessories help you to perform optical measurements without being an expert for optical measurement at all.

Illumia light sources

m-u-t has developed a series of ready-to-use light sources for the most common applications. They are designed for easy handling, equipped with SMA 905 connector and a shutter with additional external trigger control.



Light source Illumia VIS/NIR

Fiber optic cables

Quality fiber optic cables are an essential component of a spectroscopy toolbox.

m-u-t offers a whole range of fiber optic cables which are optimized for use in spectroscopy.


Sample presentation

m-u-t offers integrating spheres, cuvettes & sample holder, dip Probes and other equipment which is required to handle the light in your measuring setup.

For further information please consult the datasheet about accessories.

深圳市新世联科技有限公司

Technical specifications



Content of delivery

- Spectrometer, fully configured
- Installed entrance slit
- Calibration Protocol and data
- TRIWin 3 Pro Software
- Power supply and charger
- USB interface cable
- rugged cargo box with soft inlay

| Specifications | UV/VIS/NIR | UV/VIS | UV | VIS/NIR |
|-------------------|----------------------------------|-------------------|-------------------|-------------------|
| Spectral range | 200 - 1100 nm | 200 - 800 nm | 200 - 500 nm | 500 - 1100 nm |
| Grating | 300 l/mm | 600 l/mm | 1200 l/mm | 600 l/mm |
| Slit size (w x h) | 10 μ m x 3 mm | 40 μ m x 3 mm | 40 μ m x 3 mm | 40 μ m x 3 mm |
| Detector | back-thinned CCD, 14x 2048 pixel | | | |
| Optical design | Asymmetrical Czerny-Turner | | | |
| Connection type | SMA 905 | | | |
| Focal length | Input: 50 mm , Output: 80 mm | | | |

Performance

| | | | | |
|----------------------|---------|--------|--------|--------|
| Resolution, optical | 2.0 nm | 1.0 nm | 0.5 nm | 1.0 nm |
| Stray light | 0.1 % | 0.1 % | 0.15 % | 0.15% |
| Dark noise | 1.15 % | 1.15 % | 1.15 % | 1.15 % |
| S/N (Detector, only) | 400 : 1 | 400:1 | 400:1 | 400:1 |
| Corrected linearity | > 95 % | > 95 % | > 95 % | > 95 % |
| Wavelength Accuracy | 0.1 nm | 0.1 nm | 0.1 nm | 0.1 nm |

General

| | |
|------------------------|--|
| A/D Converter | 14 bit |
| Computer interfaces | USB, RS232 |
| Power supplies | External power supply , battery charger, rechargeable battery, NiMH 14.4 V, 1500 mAh |
| Dimensions (w x d x h) | 250 x 125 x 140 mm (10 x 5 x 6 in) |
| Weight | 2.7 kg (6 lb) |
| Onboard Processor | 32-bit RISC |
| Built-in display | Touch-screen, 256 colors |
| Triggers modes | In / out |
| Memorycard slot | Smart Media, up to 128 MB |

Part numbers

| | | |
|-----------|----------------------------|---------------|
| 30-105414 | TRISTAN 5 exUV/VIS/NIR-U-H | 200 - 1100 nm |
| 30-105413 | TRISTAN 5 exUV/VIS-U-H | 200 - 800 nm |
| 30-105410 | TRISTAN 5 exUV-U-H | 200 - 500 nm |
| 30-105997 | TRISTAN 5 VIS/NIR-U-H | 500 - 1100 nm |

